

## Lab-Tips® - LTP125 Large Knitted Polyester Swab

Cleanroom laundered double layer 100% polyester knit swab

The LTP125 large knitted polyester swab is an excellent choice for general purpose cleaning or surface sampling. Its durable flat head and rigid internal paddle provide maximum surface area for cleaning or sampling broad flat surfaces. Its cleanroom laundered knitted head provide ultra -low NVR and ion levels.



### Key Attributes

- Cleanroom laundered 100% continuous filament polyester knit head
- Long polypropylene handle with a rigid paddle head
- Thermally bonded head
- Laundered and packaged in Class 10 (ISO Class 4) cleanroom

### Benefits

- Critically low particles, fibers, ions and extractables
- High abrasion resistance while soft and non-abrasive
- Excellent chemical resistance
- No contaminating adhesives and good sorbency
- Autoclavable
- Free from Silicone, Amides and DOP

### Environmental

- Complies with RoHS and REACH requirements

### Applications

- Designed for the highest level of contamination control in critical processing applications
- Cleaning of sensitive equipment such as print heads
- Sampling surfaces for cleaning validation
- Applying lubricants and other liquids
- Removing excess materials
- Cleaning with solvents

### Physical Characteristics

<b>Head Material</b>	Knitted polyester	<b>Handle Thickness</b>	2.8mm {0.110"}
<b>Head Width</b>	12.7mm {0.50"}	<b>Handle Length</b>	100mm {3.94"}
<b>Head Thickness</b>	5mm {0.196"}	<b>Total Swab Length</b>	125mm {4.92"}
<b>Head Length</b>	25mm {0.984"}	<b>Head Bond</b>	Thermal
<b>Handle Material</b>	Polypropylene	<b>Handle Text</b>	Berkshire
<b>Handle Width</b>	5.5mm {0.216"}		

#### Other Knitted Polyester tip swabs

- LTP1465 Long handle
- LTP 70F Short handle flexible head
- LTP70R Short handle rigid head

[www.berkshire.com.sg](http://www.berkshire.com.sg)

Contact: Tel  
65.6252.4313 / Fax:  
65.6252.4312

Email: [enquiries@berkshire.com.sg](mailto:enquiries@berkshire.com.sg)

Americas	Tel 1 413 528 2602	Fax 1 413 528 2614	<a href="mailto:info@berkshire.com">info@berkshire.com</a>
Europe	Tel 44 (0) 870 757 2877	Fax 44 (0) 870 757 2878	<a href="mailto:enquiries@berkshire.uk.com">enquiries@berkshire.uk.com</a>
SE Asia	Tel 65 6252 4313	Fax 65 6252 4312	<a href="mailto:enquiries@berkshire.com.sg">enquiries@berkshire.com.sg</a>
Japan	Tel 81 3 5827 2380	Fax 81 35827 2382	<a href="mailto:master@berkshire.co.jp">master@berkshire.co.jp</a>

## Technical Data:

ATTRIBUTE		LTP125	TEST METHOD
Non-Volatile Residue	IPA mg/swab	0.075	IEST-RP.CC004.3, Sec 7.1.2
Ions (µg/swab)	Chloride	0.11	IEST-RP.CC004.3, Sec 7.2.2
	Sulfate	0.061	
	Nitrate	ND*	
	Phosphate	ND*	
	Fluoride	ND*	
	Potassium	0.010	
	Calcium	0.082	
	Sodium	0.032	
	Magnesium	0.011	

## Notes:

- Data shown are typical values
- Due to differences in test methods applied and equipment utilized by different swab manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- \* ND\* = below the detection limit of the equipment.

## Order Information:

Product	Item Number	Size	Swabs/pk	Pks/cs	Case size
Large Knitted Polyester	LTP125.5	125mm	100	5	220 x 160 X 105mm

## Other Berkshire Products




Wipers



Glove Liners



Mop Systems



Documentation Systems



Face Masks



Swabs

[www.berkshire.com.sg](http://www.berkshire.com.sg)

www.berkshire.com.sg Contact: Tel 65.6252.4313 / Fax: 65.6252.4312 Email:enquiries@berkshire.com.sg